DISPOSABLE SYRINGES AND NEEDLES


The needles are made of SUS 304 and fully coated with medical grade silicone.

Color coded hub are made of medical grade Polypropylene and transparent color to distinguish needle gauge.

Gravure printing for syringe are neat and accurate graduation serves you with precise measurement of medication.

Gasket for syringe are made of medical grade Plastic Elastomer (Latex-free) and coated with medical grade silicone for smooth extrusion.

Needle cap, barrel and plunger for syringe are made of medical grade Poly propylene compatible with any medication.





PRODUCT SPECIFICATIONS

NO	PRODUCT DESCRIPTION	SKU	QUANTITY/PACK	QUANTITY/CARTON		
Disposable Syringe With or Without Needle						
1	1 ml / cc Luer Slip (U-100)	Pack	100 pcs	800 pcs		
2	2 - 3 ml / cc Luer Slip & Luer Lock	Pack	100 pcs	800 pcs		
3	5 ml / cc Luer Slip & Luer Lock	Pack	100 pcs	800 pcs		
4	10 ml / cc Luer Slip & Luer Lock	Pack	100 pcs	800 pcs		
5	20 ml / cc Luer Slip & Luer Lock	Pack	100 pcs	800 pcs		
6	50 ml / cc Luer Lock & Luer Lock	Pack	100 pcs	800 pcs		

NO	PRODUCT DESCRIPTION	SKU	QUANTITY/PACK	QUANTITY/CARTON		
Disposable Needle						
1	Needle 19G X 1½" (40 mm)	Pack	100 pcs	3,000 pcs		
2	Needle 21 G X 1½" (40 mm)	Pack	100 pcs	3,000 pcs		
3	Needle 22 G X 1½" (40 mm)	Pack	100 pcs	3,000 pcs		
4	Needle 23 G X 1" (25 mm)	Pack	100 pcs	3,000 pcs		
5	Needle 23 G X 1¼" (30 mm)	Pack	100 pcs	3,000 pcs		
6	Needle 24 G X 1" (25 mm)	Pack	100 pcs	3,000 pcs		
7	Needle 25 G X 1" (25 mm)	Pack	100 pcs	3,000 pcs		
8	Needle 26 G X ½" (12 mm)	Pack	100 pcs	3,000 pcs		
9	Needle 26 G X 1" (25 mm)	Pack	100 pcs	3,000 pcs		

For further enquiries and details, please contact:

No. 18, Jalan Teknologi 3/5, PJU 5, Taman Sains Selangor, Kota Damansara, 47810, Petaling Jaya, Selangor Darul Ehsan, Malaysia

Tel: +603 6151 8056 Fax: +603 6143 7523 Email: info@selia-tek.com Website: www.selia-tek.com

