

INJECTION STOPPER



PRODUCT SPECIFICATIONS

Medical Grade Material

For intermittent injections through latex-free injection site

Priming volume 0.088 ml

Male Luer Lock

Luer fitting as per ISO 594-1 for compatibility with any other standard product

Minimal dead space

Individual packed in blister packaging

EtO Sterilization in accordance with ISO 11135-1

Quantity per pack: 50 pcs

Quantity per carton: 3,000 pcs

For further enquiries and details, please contact:

No. 18, Jalan Teknologi 3/5, PJU 5, Taman Sains Selangor, Kota Damansara,
47810, Petaling Jaya, Selangor Darul Ehsan, Malaysia

Tel: +603 6151 8056 **Fax:** +603 6143 7523

Email: info@selia-tek.com **Website:** www.selia-tek.com