















SELIA-TEK
moulding the future

VACTUBE

BLOOD COLLECTION TUBE



VACTUBE ORDER OF DRAW

ADDITIVE/CAP COLOUR	RECOMMENDED INVERTIONS TIME	TYPE OF TEST
No Additive	0	To use as discard tube.
		
3.2% Sodium Citrate	3 - 4	For routine coagulation tests or selected platelet function assays.
		
Clot Activator	5 - 8	For blood chemistry assays, routine blood donor screening infectious disease testing.
		
Clot Activator with Gel	5 - 8	Contain gel for serum separation; for blood chemistry assays, hormone studies, routine blood donor screening, infectious disease testing.
		
Lithium Heparin	5 - 8	For routine blood chemistry tests, cardiac enzymes tests, bone profile testing.
		
Sodium Heparin	5 - 8	For routine blood chemistry tests, cardiac system enzymes tests.
		
Lithium Heparin with Gel	5 - 8	Contain gel for serum separation for routine blood chemistry tests, cardiac enzymes tests, bone profile testing.
		
EDTA K2	5 - 8	For whole blood haematology study, routine immunohematology tests, blood donor screening, cross match profile.
		
EDTA K3 with Gel	5 - 8	Contain gel for plasma separation; for whole blood haematology study, routine immunohematology test, blood donor screening.
		
Sodium Fluoride	5 - 8	For glucose and lactate testing.
		
ESR	8 - 10	For Erythrocyte Sedimentation Rate (ESR) test.
		

**Product colour may slightly vary due to photographic lighting sources

VACTUBE PLAIN

The tube contains no additive.



REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
10040B	13x75	4.0ml	No Additive	18 months	100x10

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
11060B	13x100	6.0ml	No Additive	18 months	100x10

VACTUBE SODIUM CITRATE

The tube contains 3.2% Sodium Citrate as additive. The tube is made up of TWO (2) layers: Polypropylene (PP) for inner tube and Polyethylene Terephthalate (PET) as the outer tube.

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
200401	13x75	1.8ml	Sodium Citrate	18 months	100x10

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
210601	13x75	2.7ml	Sodium Citrate	18 months	100x10



**Product colour may slightly vary due to photographic lighting sources

VACTUBE CLOT ACTIVATOR

The tube contains Silica Particles as additive.



REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
200401	13x75	4.0ml	Silica Particles	18 months	100x10

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
210601	13x100	6.0ml	Silica Particles	18 months	100x10

VACTUBE CLOT ACTIVATOR WITH GEL

The tube contains Silica Particles as additive and get that act as serum separator.

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
201352	13x75	3.5ml	Silica Particles	18 months	100x10

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
211502	13x100	5.0ml	Silica Particles	18 months	100x10

**Product colour may slightly vary due to photographic lighting sources



VACTUBE SODIUM HEPARIN

The tube contains Sodium Heparin as additive.



REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
200402	13x75	4.0ml	Sodium Heparin	18 months	100x10

VACTUBE LITHIUM HEPARIN

The tube contains Lithium Heparin as additive.

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
200401	13x75	4.0ml	Lithium Heparin	18 months	100x10

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
210601	13x100	6.0ml	Lithium Heparin	18 months	100x10



**Product colour may slightly vary due to photographic lighting sources

VACTUBE LITHIUM HEPARIN WITH GEL

The tube contains Lithium Heparin as additive and gel that act as plasma separator.



REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
601306	13x75	3.0ml	Lithium Heparin	18 months	100x10

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
601356	13x75	3.5ml	Lithium Heparin	18 months	100x10

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
611506	13x100	5.0ml	Lithium Heparin	18 months	100x10

VACTUBE EDTA K2

The tube contains K2 EDTA as additive.

EDTA - Ethylenediaminetetraacetic acid

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
400203	13x75	2.0ml	EDTA K2	18 months	100x10

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
400303	13x75	3.0ml	EDTA K2	18 months	100x10

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
400403	13x75	4.0ml	EDTA K2	18 months	100x10

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
410603	13x100	6.0ml	EDTA K2	18 months	100x10



**Product colour may slightly vary due to photographic lighting sources

VACTUBE EDTA K3 WITH GEL

The tube contains EDTA K3 as additive and gel that act as plasma separator.

EDTA - Ethylenediaminetetraacetic acid



REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
301354	13x75	3.5ml	EDTA K3	18 months	100x10

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
311504	13x100	5.0ml	EDTA K3	18 months	100x10

VACTUBE SODIUM FLUORIDE

The tube contains Sodium Fluoride and EDTA K3 as additive.

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
800208	13x75	2.0ml	Sodium Fluoride & EDTA K3	18 months	100x10

REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
800408	13x75	4.0ml	Sodium Fluoride & EDTA K3	18 months	100x10

**Product colour may slightly vary due to photographic lighting sources



VACTUBE ESR

The tube contains 3.8% Sodium Citrate as additive for Erythrocyte Sedimentation Rate (ESR) test.



REF. NUMBER	TUBE SIZE	DRAW VOLUME	ADDITIVE	SHELF LIFE	PACKAGING SIZE (CARTON)
A00169	13x75	1.6ml	3.8% Sodium Citrate	18 months	100x10

**Product colour may slightly vary due to photographic lighting sources

NOTES

Awards & Achievements



MINOLTA

Best Delivery Award by Minolta - 2001



ISO MS9001:2000 certification - 2003



Golden Bull Awards Top 100 SMIs - 2003



Enterprise 50 Awards (27th placed) - 2003



TNB Mentor and Mentee Program - 2004



Appointed as TNB Preferred Vendor - 2007



Appointed by MPC as "Model Company for Development of Excellent Organisation" - 2007



ISO 13485:2016 certification - 2007



5S Certification by MPC - 2008



ISO 9001:2015 Certification - 2009



Selected to participate in TQM Model Company Conducted by APO/MPC - 2009



Selected under MATRADE Bumiputra Entrepreneur Programme - 2010



Successfully obtaining 3 stars in SME COMPETITIVENESS RATING FOR ENHANCEMENT (SCORE) by SME Corporation - 2010



Selected to participate in Coaching & Mentoring Programme By Technology Park Malaysia Corporation under MITI - 2011



Enterprise 50 Awards (49th placed) - 2012



Successfully obtaining 4 stars in SME COMPETITIVENESS RATING FOR ENHANCEMENT (SCORE) by SME Corporation - 2012



Successfully obtaining 4 stars in SME COMPETITIVENESS RATING FOR ENHANCEMENT (SCORE) by SME Corporation - 2013



Enterprise 50 Awards (19th placed) - 2013



SME 100 Award : Fast Moving Companies - 2014



Selected to Participate 1-Innocert - 2014



MS ISO/IEC 17021 - 2006 QS 091220912 CB 12



015



Medical Products| SeliaTek

For further enquiries and details, please contact:



Selia-Tek Holdings Sdn Bhd (141025-X)

Lot.18, Jalan Teknologi 3/5, PJU 5, Taman Sains Selangor,
47810 Kota Damansara, Selangor Darul Ehsan, Malaysia

Tel: +603 6151 8056 - 9

Fax: +603 6150 1771

Email: info@seliatek.com

Website: www.seliatek.com

- **Selia-Tek Industries Sdn Bhd** (199531-V)
- **Selia-Tek Medical Sdn Bhd** (567990-H)

- **Selia-Tek Ventures Sdn Bhd** (740907)
- **Selia Bio Medic Sdn Bhd** (861451-A)